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Dr Ana G Cabado
Food Safety Division
ANFACO-CECOPECA

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





11.1.2011

EN

Official Journal of the European Union

L 6/3

COMMISSION REGULATION (EU) No 15/2011

of 10 January 2011

amending Regulation (EC) No 2074/2005 as regards recognised testing methods for detecting marine biotoxins in live bivalve molluscs

(Text with EEA relevance)


EU-RL-MB
EUROPEAN UNION REFERENCE
LABORATORY FOR MARINE
BIOTOXINS

**EU-Harmonised Standard Operating
Procedure for determination of
Lipophilic marine biotoxins in molluscs
by LC-MS/MS**

Version 5, January
2015

Coordination:
European Union Reference Laboratory for Marine Biotoxins (EU-RL-MB)
Agencia Española de Consumo, Seguridad Alimentaria y Nutrición (AECOSAN)
CITEXVI (Ciudad Tecnológica de Vigo)
Fonte das Abelieiras, s/n, Campus Universitario de Vigo, As Lagoas-Marcoense
36310 Vigo, Spain
e-mail: eurimb@csssi.es
<http://aesan.msssi.gob.es/en/CRLMB/web/home.shtml>

Plan de la page Web

European Union Reference Laboratory for Marine Biotoxins

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Other Standard Operating Procedures

METHOD FOR QUANTITATIVE DETERMINATION OF MARINE LIPOPHILIC TOXINS IN SHELLFISH USING LC-MS/MS

Agosto 2011

Standard Operating Procedures method for the quantitative determination of marine lipophilic toxins in shellfish using LC-MS/MS.

[Download Document](#)

H.J. van den Top, A. Gerssen, P. McCarron, H.P. van Egmond. Quantitative determination of marine lipophilic toxins in mussels, oysters and cockles using liquid chromatography-mass spectrometry: inter-laboratory validation study. Food Additives & Contaminants: Part A, 2011, Vol. 28, Iss. 12.

METHOD FOR DETECTION OF OKADAIC ACID TOXINS GROUP IN EDIBLE MOLLUSCS Rev 4

October 2010

Standard Operating Procedures method for the quantitative determination of OA and other carboxylic toxins of the OA group including DTX1, DTX2, DTX3 by a colorimetric phosphatase inhibition assay.

The OkaTest method fulfils the requirements of Chapter III A (4) a, b, and c of Annex III to Regulation (EC) 2074/2005 and can be used as supplementary method to EU-RL LC/MS/MS method only for the detection of Okadaic Acid group toxins OA, DTX1, DTX2, DTX3 including their esters.

The use of OkaTest method is not an alternative to the EU-RL LC/MS/MS method.

LC-MS/MS conditions

Column	Waters Acquity UPLC BEH C18 1.7 μ m 2.1 x 100 mm	
Flow	0.40 mL / min	
Injection volume	2 μ l	
Column temperature	45 $^{\circ}$ C	
Maximum pressure	15000 psi	
Time (minutes)	% Mobile phase A Ammonia 6.7 mM in water	% Mobile phase B Ammonia 6.7 mM in acetonitrile
initial	70	30
0.50	70	30
6.00	10	90
6.50	0	100
9.50	0	100
9.60	70	30
10.00	70	30

MS-MS fragmentation

Compound	Parent ion (Da)	Fragment ion (Da)	ESI	Cone voltage (V)	Collision voltage (eV)	Dwell (s)
AO-DTX2	803.50	113.00	-	86	56	0.007
AO-DTX2*	803.50	255.20	-	86	46	0.007
DTX1*	817.50	255.20	-	86	48	0.007
DTX1	817.50	113.1	-	86	43	0.007
DTX1	817.50	113.1	-	86	43	0.007

UPLC (Acquity H-class) acoplado a un detector MS/MS detector (XEVO TQS) y a un Acquity UPLC BEH C18 1.7 μ m 2.1 \times 100 mm column (Waters)

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Toxinas ASP toxins: Domoic acid (DA)

CHAPTER II

AMNESIC SHELLFISH POISON (ASP) DETECTION METHOD

The total content of amnesic shellfish poison (ASP) of edible parts of molluscs (the entire body or any part edible separately) must be detected using the high-performance liquid chromatography (HPLC) method or any other recognised method.

If the results are challenged, the reference method shall be the HPLC method.

22.12.2005

EN

Official Journal of the European Union

L 338/27

COMMISSION REGULATION (EC) No 2074/2005

of 5 December 2005

laying down implementing measures for certain products under Regulation (EC) No 853/2004 of the European Parliament and of the Council and for the organisation of official controls under Regulation (EC) No 854/2004 of the European Parliament and of the Council and Regulation (EC) No 882/2004 of the European Parliament and of the Council, derogating from Regulation (EC) No 852/2004 of the European Parliament and of the Council and amending Regulations (EC) No 853/2004 and (EC) No 854/2004

HPLC-UV

EUROPEAN UNION REFERENCE LABORATORY
FOR MARINE BIOTOXINS

**EU-Harmonised Standard Operating
Procedure for determination of
domoic acid in shellfish and finfish by
RP-HPLC using UV detection**

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**EUROPEAN UNION REFERENCE
LABORATORY FOR MARINE BIOTOXINS**

**EURLMB Standard Operating Procedure for
PSP toxins by Mouse Bioassay**

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JRC TECHNICAL REPORTS

Algal bloom and its economic impact


Figure 4: Total Number of Harmful Algae Events reported globally by single countries during the period 1980-2015 to the Harmful Algal Events Dataset (HAEDAT)

Isabella Sanseverino, Diana Conduto, Luca Pozzoli, Srdan Dobricic and Teresa Lettieri

2016

Number of Harmful Algae Events by Syndrome




Figure 5: Total Number of Harmful Algae Events by syndrome reported globally during the period 1980-2015 to the Harmful Algal Events Dataset (HAEDAT)

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


Toxicon 55 (2010) 235–243

Contents lists available at ScienceDirect

Toxiconjournal homepage: www.elsevier.com/locate/toxicon

Decrease of marine toxin content in bivalves by industrial processes

Antonio Reboreda^a, Jorge Lago^a, María-José Chapela^a, Juan M. Vieites^a, Luis M. Botana^b, Amparo Alfonso^b, Ana G. Cabado^{a,*}

^aMicrobiology and Toxins Area, ANFACO-CECOPESCA, Campus Univ. 16, 36310 Vigo PO, Spain

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
Fig. 4. Effect of thermal processing or/and freezing on PSP toxin levels in clams ($n=4$). Values obtained were $700 \pm 60 \mu\text{g eq STX/kg}$ for fresh clams and $873.3 \pm 17.6 \mu\text{g eq STX/kg}$ for frozen clams. Thermal processing decreased toxicity below the levels of detection.


Fig. 5. Effect of thermal processing or/and freezing in PSP toxin levels on cockles ($n=4$). Control PSP concentration was $673.3 \pm 3.3 \mu\text{g eq STX/kg}$. Thawed cockles showed a significant reduction on PSP concentration $366.70 \pm 16.70 \mu\text{g eq STX/kg}$. Thermal processing reduced toxicity below the detection limit.

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An aerial photograph of a large reservoir. In the background, a suspension bridge spans across the water. The water is a deep blue-green color, and the surrounding hills are covered in dense green forest. In the foreground, there is a small island with a few buildings and trees, connected to the mainland by a small bridge. The text "Gracias por su atención" is overlaid in the center of the image.

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